



Discover a tube whose quality speaks for itself

Velocity™ Single Buccal Tubes


It's the newest design of an extremely low-profile buccal tube, and patient comfort and clinical effectiveness are the results. Funneled entrance, occlusal/gingival positioning guides, distal notch for additional ligation. 4.3mm.

Blaze™ Curing Light

- **Added Feature:** Caries detection light
- Colorful, ergonomic handle design
- Color choices: Grass Green, Pink, Purple, Blue and Black
- A full charge will provide more than 200 applications
- Three working modes: full power, ramping-up and pulse
- Variable time settings: from 5 seconds to 40 seconds
- Long stand-by time, as long as six months
- Wave length: 420nm - 480nm
- Light intensity: 1,600 mW/cm² to 1,800 mW/cm²
- Specify Black, Purple, Blue, Green or Pink Accents
- 1 year warranty

\$199 ea.

SAVE 50%



Tweezer notch and flat buccal surface for easy handling and seating

Generously sized and contoured pad with indent for error-free placement and bond reliability

Torque-in base for accurate prescription expression

ORDER TODAY!
Items must ship by August 31st

Call (800) 854-2896

Special Limited-Time Pricing \$4.75 ea.

Buy 100 get 20 Free ... Under \$4.00 ea.

Buy 160 get 40 Free ... SAVE 20%

Buy 300 get 100 Free ... SAVE 25%

Add Band Assemblies for \$1.99 ea.

LANCER
ORTHODONTICS



50th
Anniversary

EST. 1967