


orthotown[®]

A division of Farran Media

real orthodontics for real orthodontists

July/August 2018 • Vol. 11, Issue 6



TOWNIE CHOICE AWARDS PRODUCT SHOWCASE PREVIEW

p. 42



ANCHORS AWEIGH
Dr. Saritha Chary-Reddy and Karthik Reddy share a microscrew-assisted Class II treatment case **p. 36**




**ADULT PATIENTS
IN CLEAR ALIGNERS**
Dr. Darren Forcier discusses aesthetics-driven treatment **p. 28**



LANCER ORTHODONTICS

Praxis



Orthodontists choose Praxis when treatment requires an exceptionally well-designed bracket that promises—and delivers—the highest level of precision. Praxis brackets have become a favorite among orthodontists worldwide.

The Praxis mini-twin combines robotic metal injection-molded technology, superior design and pleasing aesthetics to deliver a robust bracket that is inherently stronger and more dimensionally accurate than similar brackets. The bracket is loaded with features that patients and doctors love, including precise, diamond-milled slots. In addition, the torque-in-base allows for fewer wire bends and repositioned brackets. A Scribe line runs through the bracket and the pad to facilitate fast and accurate placement.

Lancer's Praxis brackets turn theory into reality, combining art, science and skill to create an easy-to-use, beautifully functional bracket that consistently produces excellent treatment results. Manufactured under a stringent, FDA- and ISO-regulated quality management system that helps ensure only the highest-quality parts are delivered to its customers, the Praxis bracket stands out among the best in the industry in terms of design, performance and overall quality. The brackets are available in Lancer's patented Redi-Pak,

which allows clinicians to order and store their brackets in a sealed package until patients are ready for bonding. ■

HIGHLIGHTS

- **Precision-milled slots** guarantee precise dimensions throughout the wall of the wire slot.
- **Large under-tie-wing area** for easy ligation.
- **Optimum low-profile design** for optimum patient comfort.

For more information, visit lancerortho.com or call 800-854-2896.

